



# BeMicron® 4AII

Classic cleanroom coverall with a separate hood. Two models to choose from: stand-up collar and knit collar.



(· 🗸 ) 30 sec

**ESD** 

Grades A,B,C, D

ISO 3 → 7

The BeMicron 4All coverall offers a modern gowning system in a classic design. It can be used in virtually all types of clean room.

Made from a fabric carefully selected for its quality, performance and lasting reliability, it offers excellent mechanical strength and protects against contamination. It has a perfectly-smooth technical filtering wall, at both front and back, and constitutes a veritable barrier against human contamination. Thanks to the BeMicron® folding system and donning procedure, the operator's movements follow each other automatically without need to touch the outside of the coverall and without errors in procedure.

## The doning procedure



Out of the bag



White = left Red = right



Easy gowning system



**Easely** adjustable



Zipper quickly

30 sec







### BeMicron® 4All



#### Classical officer neck: closure by 3 stainless steel studs aligned with collar bone

In order to allow improved and comfortable gowning (passing of the head), studs have been added.



# Knitted neck: easy to wear and adjusts automatically

The knitted collar allows easy gowning and automatic adjustment around the neck.



## **Durable and effective** cuff

The cuff or stud finish allows the sleeves to fit perfectly into gloves. The cuff is made of antistatic polyester and is resistant to various sterilization processes. The studs are made of stainless steel.



# Protective clothing with optimum comfort

The elastic in the back ensures the BeMicron adjusts perfectly to the body to prevent any «blow-up effect» caused by movements in the body.



## Gowning without touching the coverall

The BeMicron is unique in the fact that the outside of the garment is protected from any risk of contact contamination (hand and surface). The base of the zip is already engaged and in place reducing contact time.. Simply use the grip to close the zip.

